

Product Datasheet

PTH antibody (orb76758)

Description	Mouse monoclonal antibody to PTH
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, RIA
Immunogen	Synthetic PTH (aa1-34) conjugated to carrier protein
Target	Parathyroid Hormone (PTH) N-terminal
Preservatives	Purified antibody (from supernatant) containing PBS + 0.1% sodium azide
Concentration	2mg/ml
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
Uniprot ID	P01270
Entrez	5741
Dilution Range	Optimal antibody dilution should be determined by titration
Expiration Date	12 months from date of receipt.