



Product Datasheet

CKR-5 (phospho-Ser336) antibody (orb767585)

Description	Rabbit polyclonal antibody to CKR-5 (phospho-Ser336)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human CCR5 around the phosphorylation site of Ser336. AA range:302-351
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	40
Expiration Date	12 months from date of receipt.



Western
Blot
analysis of
Jurkat cells
US...



Western
blot
analysis of
lysates from
Ju...