



Product Datasheet

CB1 antibody (orb767598)


Description	Rabbit polyclonal antibody to CB1
Species/Host	Rabbit
Reactivity	Human, Monkey, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human CNR1. AA range:151-200
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	53
Expiration Date	12 months from date of receipt.



Western Blot analysis of HT-29 cells using CB1 antibody.



Immunofluorescence analysis of LOVO cells using CB1 antibody.



Immunohistochemistry analysis of paraffin-embedded tissue using CB1 antibody.