


Product Datasheet


COL11A1 antibody (orb767627)

Description	Rabbit polyclonal antibody to COL11A1
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human Collagen XI alpha1. AA range:581-630
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	181
Expiration Date	12 months from date of receipt.




Immunofluorescence analysis of A549. 1,

...



Western Blot analysis of KB cells using ...



Western blot analysis of lysates from K5...