



Product Datasheet

ADORA2A antibody (orb767652)


Description	Rabbit polyclonal antibody to ADORA2A
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human ADORA2A. AA range:120-169
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	37
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of COS7 cell...



Immunohistochemistry analysis of paraffi...



Western blot analysis of ADORA2A Antibod...