


Product Datasheet

PDE7A antibody (orb76771)

Description	Rabbit monoclonal antibody to PDE7A
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, IP, WB
Immunogen	Synthetic peptide, corresponding to amino acids at the C-terminal of Human PDE7A
Target	PDE7A
Preservatives	Diluted serum
Note	For research use only
Clonality	Monoclonal
Uniprot ID	O08593 , P70453 , Q13946
Entrez	5150 , 81744 , 18583
Dilution Range	Optimal antibody dilution should be determined by titration, however as a guideline try; WB: 1/2000. Predicted molecular weight: 56 kDa
Expiration Date	12 months from date of receipt.



Western blot analysis of untransfected u...