



Product Datasheet

Casein Kinase Ialpha (phospho-Tyr321) antibody (orb767731)


Description	Rabbit polyclonal antibody to Casein Kinase Ialpha (phospho-Tyr321)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human Casein Kinase I alpha around the phosphorylation site of Tyr321. AA range:287-336
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	35
Expiration Date	12 months from date of receipt.



Western Blot analysis of VEC cells using...



Enzyme-Linked Immunosorbent Assay (Phosp...



Western blot analysis of lysates from HU...