



Product Datasheet

Casein Kinase IEpsilon antibody (orb767734)


Description	Rabbit polyclonal antibody to Casein Kinase IEpsilon
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human CKI-epsilon. AA range:276-325
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	47
Expiration Date	12 months from date of receipt.



Western Blot analysis of LOVO cells using...



Immunofluorescence analysis of HeLa cells...



Immunohistochemistry analysis of paraffin...