



Product Datasheet

Six5 antibody (orb767752)


Description	Rabbit polyclonal antibody to Six5
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human SIX5. AA range:201-250
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	75
Expiration Date	12 months from date of receipt.



Western Blot analysis of K562 cells using Six5 antibody.



Western blot analysis of lysates from K562 cells using Six5 antibody.



Immunohistochemical analysis of paraffin-embedded tissue sections using Six5 antibody.