


Product Datasheet

PP2A C antibody (orb76780)

Description	Mouse monoclonal antibody to PP2A C
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	A KLH-coupled peptide with the sequence: KLH-C- $\ddot{A}$ YA-RRTPDYFL-OMe. The c-terminal sequence (as well as the carboxy-methylation) is absolutely conserved from human, mouse to yeast
Target	PP2A alpha (methyl L309)
Preservatives	Purified antibody (from supernatant) containing PBS + 0.1% sodium azide
Concentration	1mg/ml
Note	For research use only
Application notes	X63-Ag8.653
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	P63330, P63331, P67775
Entrez	5515, 24672, 19052
Dilution Range	Optimal antibody dilution should be determined by titration, however as a guideline try; WB: 1:500. Detects a band of approximately 36 kDa (predicted molecular weight: 36 kDa)
Expiration Date	12 months from date of receipt.



Western blot analysis of native lysate
u...