



Product Datasheet

Pragmin (phospho-Tyr413) antibody (orb767817)


Description	Rabbit polyclonal antibody to Pragmin (phospho-Tyr413)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	Synthesized phospho-peptide around the phosphorylation site of human Pragmin (phospho-Tyr413)
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	175
Expiration Date	12 months from date of receipt.



Western Blot analysis of A549 cells using...



Enzyme-Linked Immunosorbent Assay (Phosp...



Western blot analysis of lysate from A54...