



Product Datasheet

AKAP 14 antibody (orb767824)


Description	Rabbit polyclonal antibody to AKAP 14
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human AKAP14. AA range:1-50
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	23
Expiration Date	12 months from date of receipt.



Western Blot analysis of Jurkat cells us...



Immunofluorescence analysis of HeLa cell...



Western blot analysis of lysates from Ju...