



Product Datasheet

ACE1 antibody (orb767848)


Description	Rabbit polyclonal antibody to ACE1
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human ACE1. AA range:891-940
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	165
Expiration Date	12 months from date of receipt.



Western Blot analysis of HuvEc cells using ACE1 antibody.



Immunohistochemistry analysis of ACE1 antibody in tissue sections.



Western blot analysis of lysates from mouse cells using ACE1 antibody.