

Product Datasheet

Progesterone Receptor C262 antibody (orb76787)

Description	Mouse monoclonal antibody to Progesterone Receptor C262
Species/Host	Mouse
Reactivity	Human, Mouse, Rabbit, Rat
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	Synthetic peptide to C-terminal tail of human progesterone receptor (PR). (aa 919-933 ILAGMVKPLLHKK)
Target	Progesterone Receptor
Preservatives	Purified antibody containing PBS 0.1% sodium azide
Concentration	1mg/ml
Note	For research use only
Application notes	FoxNY (mouse myeloma from Balb/c origin)
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	Q63449 , P06401 , Q00175 , P06186
Entrez	18667 , 100009094 , 5241 , 25154
Dilution Range	Optimal antibody dilution should be determined by titration, however as a guideline use at 5-10ug/ml
Expiration Date	12 months from date of receipt.