



Product Datasheet

Photomedin-1 antibody (orb767871)


Description	Rabbit polyclonal antibody to Photomedin-1
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human OLFML2A. AA range:311-360
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	73
Expiration Date	12 months from date of receipt.



Western Blot analysis of RAW264.7 cells ...



Immunofluorescence analysis of A549 cell...



Western blot analysis of lysates from RA...