



Product Datasheet

TdT antibody (orb767889)


Description	Rabbit polyclonal antibody to TdT
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human DNTT. AA range:381-430
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	60
Expiration Date	12 months from date of receipt.



Western Blot analysis of HepG2 cells using TdT antibody.



Immunohistochemical analysis of paraffin-embedded tissue using TdT antibody.



Immunohistochemical analysis of paraffin-embedded tissue using TdT antibody.