

Product Datasheet

Progesterone Receptor N559 antibody (orb76790)

Description	Mouse monoclonal antibody to Progesterone Receptor N559
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	Synthetic peptide to human progesterone receptor (PR). Amino-terminal domain aa. 551-564 (ASQSPQTSFESLPQ)
Target	Progesterone Receptor
Preservatives	Purified antibody containing PBS + 0.1% sodium azide
Concentration	1mg/ml
Note	For research use only
Application notes	FoxNY (mouse myeloma from Balb/c origin)
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	Q00175
Entrez	18667
Dilution Range	Optimal antibody dilution should be determined by titration, however as a guideline use at 5-10ug/ml
Expiration Date	12 months from date of receipt.