



Product Datasheet

D1DR antibody (orb767901)


Description	Rabbit polyclonal antibody to D1DR
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human DRD1. AA range:135-184
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	50
Expiration Date	12 months from date of receipt.



Western Blot analysis of HepG2 cells using...



Immunohistochemical analysis of paraffin-embedded tissue sections...



Immunofluorescence analysis of MCF7 cells...