



Product Datasheet

AT1 antibody (orb767907)


Description	Rabbit polyclonal antibody to AT1
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human AGTR1. AA range:101-150
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	41
Expiration Date	12 months from date of receipt.



Western Blot analysis of HEPG2 cells using AT1 antibody.



Immunofluorescence analysis of MCF7 cells using AT1 antibody.



Western blot analysis of lysates from K562 cells using AT1 antibody.