



Product Datasheet

Ephrin-B1 (phospho-Tyr317) antibody (orb767920)


Description	Rabbit polyclonal antibody to Ephrin-B1 (phospho-Tyr317)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human EFNB1 around the phosphorylation site of Tyr317. AA range:283-332
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	30
Expiration Date	12 months from date of receipt.



Western Blot analysis of mouse cells using...



Enzyme-Linked Immunosorbent Assay (Phosp...



Western blot analysis of lysates from mouse cells...