



Product Datasheet

EphA2/3/4 (phospho-Tyr588/596) antibody (orb767947)


Description	Rabbit polyclonal antibody to EphA2/3/4 (phospho-Tyr588/596)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human EPHA2/3 around the phosphorylation site of Tyr588/596. AA range:556-605
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	130
Expiration Date	12 months from date of receipt.



Western Blot analysis of HepG2 cells using...



Enzyme-Linked Immunosorbent Assay (Phospho...)



Immunofluorescence analysis of HeLa cells...