


## Product Datasheet

### EphA2/5 (phospho-Tyr594) antibody (orb767948)

|                            |                                                                                               |
|----------------------------|-----------------------------------------------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to EphA2/5 (phospho-Tyr594)                                        |
| <b>Species/Host</b>        | Rabbit                                                                                        |
| <b>Reactivity</b>          | Human, Mouse                                                                                  |
| <b>Conjugation</b>         | Unconjugated                                                                                  |
| <b>Tested Applications</b> | ELISA, WB                                                                                     |
| <b>Immunogen</b>           | Synthesized phospho-peptide around the phosphorylation site of human EphA2/5 (phospho Tyr594) |
| <b>Preservatives</b>       | PBS with 0.02% sodium azide and 50% glycerol pH 7.4.                                          |
| <b>Concentration</b>       | 1 mg/ml                                                                                       |
| <b>Storage</b>             | Store at -20°C. Avoid repeated freeze-thaw cycles.                                            |
| <b>Note</b>                | For research use only                                                                         |
| <b>Isotype</b>             | IgG                                                                                           |
| <b>Clonality</b>           | Polyclonal                                                                                    |
| <b>MW</b>                  | 110                                                                                           |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                               |



Western  
Blot  
analysis of  
JK cells  
using ...