


Product Datasheet

RanGAP antibody (orb76796)

Description	Rabbit monoclonal antibody to RanGAP
Species/Host	Rabbit
Reactivity	Human, Mouse, Xenopus
Conjugation	Unconjugated
Tested Applications	ICC, IHC-P, WB
Immunogen	Synthetic peptide mixture: CHWSDMFTGRLRPEI (amino acids 82-96) and CPSPEKLVRMGPRRSA (amino acids 434-448) from <i>Xenopus laevis</i>
Target	RanGAP1
Preservatives	Antiserum containing PBS 0.1% sodium azide
Note	For research use only
Clonality	Monoclonal
Uniprot ID	P46061 , O13066 , P46060
Entrez	19387 , 5905
Dilution Range	Optimal antibody dilution should be determined by titration however as a guideline try; WB: 1/500 - 1/1000. Detects a band of approximately 60 kDa
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...