



Product Datasheet

HuB/D antibody (orb767962)


Description	Rabbit polyclonal antibody to HuB/D
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human ELAV2/4. AA range:61-110
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	40
Expiration Date	12 months from date of receipt.



Western Blot
analysis of HeLa cells
usin...



Immunofluorescence
analysis of HepG2
cel...



Western blot
analysis of lysates
from Ju...