



Product Datasheet

Elf-1 antibody (orb767963)


Description	Rabbit polyclonal antibody to Elf-1
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human ELF1. AA range:537-586
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	67
Expiration Date	12 months from date of receipt.



Western Blot
analysis of 293 cells
using...



Immunofluorescence
analysis of A549
cell...



Western blot
analysis of lysates
from 29...