



Product Datasheet

Neu (phospho-Tyr1248) antibody (orb768016)


Description	Rabbit polyclonal antibody to Neu (phospho-Tyr1248)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human HER2 around the phosphorylation site of Tyr1248. AA range:1206-1255
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	180
Expiration Date	12 months from date of receipt.



Western Blot analysis of 293 cells using...



Western Blot analysis of HepG2 + PMA cel...



Immunohistochemistry analysis of paraffi...