

Product Datasheet

Retinal Santigen antibody (orb76802)

Description Mouse monoclonal antibody to Retinal Santigen

Species/Host Mouse

Reactivity Bovine, Human, Porcine, Rat

Conjugation Unconjugated

Tested Applications ICC, IHC-Fr, WB

Immunogen Porcine S-antigen (arrestin)

Target Retinal S antigen

Preservatives Purified antibody (from supernatant) containing PBS 0.1% sodium azide

Concentration 0.6mg/ml

Note For research use only

Application notes Balb/c Ag 8.653 myeloma cells

Isotype IgG1


Clonality Monoclonal

Uniprot ID **P10523**

Entrez **6295**

Dilution Range Optimal antibody dilution should be determined by titration, however as a guideline try at 1:100 for IHC (on bovine retina)

Expiration Date 12 months from date of receipt.



Immunohistochemical analysis of formalin...