



Product Datasheet

ErbB-3 (phospho-Tyr1289) antibody (orb768028)


Description	Rabbit polyclonal antibody to ErbB-3 (phospho-Tyr1289)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human HER3 around the phosphorylation site of Tyr1289. AA range:1256-1305
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	185
Expiration Date	12 months from date of receipt.



Western Blot analysis of HeLa cells using anti-ErbB-3 (phospho-Tyr1289) antibody.



Enzyme-Linked Immunosorbent Assay (Phospho-ELISA) showing high specificity for the phosphopeptide.



Immunofluorescence analysis of HeLa cells using anti-ErbB-3 (phospho-Tyr1289) antibody.