



Product Datasheet

Akt (phospho-Ser124) antibody (orb768034)


Description	Rabbit polyclonal antibody to Akt (phospho-Ser124)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human Akt around the phosphorylation site of Ser124. AA range:90-139
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	55
Expiration Date	12 months from date of receipt.



Western Blot analysis of mouse cells using Akt (phospho-Ser124) antibody (orb768034).



Immunohistochemistry analysis of paraffin-embedded tissue sections using Akt (phospho-Ser124) antibody (orb768034).



Western blot analysis of lysates from mouse cells using Akt (phospho-Ser124) antibody (orb768034).