



Product Datasheet

Cleaved-Factor XII HC (I20) antibody (orb768065)

Description	Rabbit polyclonal antibody to Cleaved-Factor XII HC (I20)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human FA12. AA range:1-50
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.



Western blot analysis of FA12 (heavy cha...



Immunohistochemistry analysis of paraff...