

Product Datasheet

PYK2 (phospho-Tyr402) antibody (orb768074)

Description Rabbit polyclonal antibody to PYK2 (phospho-Tyr402)

Species/Host Rabbit

Reactivity Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications ELISA, IHC-P

Immunogen The antiserum was produced against synthesized peptide derived from human PYK2 around the phosphorylation site of Tyr402. AA range:369-418

Preservatives PBS with 0.02% sodium azide and 50% glycerol pH 7.4.

Concentration 1 mg/ml


Storage Store at -20°C. Avoid repeated freeze-thaw cycles.

Note For research use only


Isotype IgG

Clonality Polyclonal

Expiration Date 12 months from date of receipt.



Immunohistochemistry analysis of paraffi...



Western blot analysis of PYK2 (Phospho-T...