



Product Datasheet

IF3 antibody (orb768080)


Description	Rabbit polyclonal antibody to IF3
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human MTIF3. AA range:69-118
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	30
Expiration Date	12 months from date of receipt.



Western Blot analysis of HuvEC cells using IF3 antibody.



Western blot analysis of the lysates from MTIF3-expressing cells.



Immunohistochemical analysis of paraffin-embedded tissue sections using IF3 antibody.