



Product Datasheet

Olfactory receptor 5M3 antibody (orb768088)


Description	Rabbit polyclonal antibody to Olfactory receptor 5M3
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human OR5M3. AA range:53-102
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	35
Expiration Date	12 months from date of receipt.



Western blot analysis of the lysates fro...



Western Blot analysis of HEPG2 cells usi...



Immunofluorescence analysis of MCF7 cell...