



Product Datasheet

Olfactory receptor 4D6 antibody (orb768096)


Description	Rabbit polyclonal antibody to Olfactory receptor 4D6
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human OR4D6. AA range:260-309
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	36
Expiration Date	12 months from date of receipt.



Western
Blot
analysis of
LOVO cells
usin...



Western
blot
analysis of
lysates from
LO...



Western
blot
analysis of
the lysates
fro...