

Product Datasheet

Intestinal Cell Kinase (phospho-Tyr159) antibody (orb768137)

Description Rabbit polyclonal antibody to Intestinal Cell Kinase (phospho-Tyr159)

Species/Host Rabbit

Reactivity Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications ELISA, WB

Immunogen The antiserum was produced against synthesized peptide derived from human ICK around the phosphorylation site of Tyr159. AA range:125-174

Preservatives PBS with 0.02% sodium azide and 50% glycerol pH 7.4.

Concentration 1 mg/ml

Storage Store at -20°C. Avoid repeated freeze-thaw cycles.


Note For research use only

Isotype IgG


Clonality Polyclonal

MW 71


Expiration Date 12 months from date of receipt.



Western Blot analysis of 3T3 cells using...



Western blot analysis of lysates from NI...



Immunohistochemical analysis of paraffin...