



Product Datasheet

Shank 2 antibody (orb768152)


Description	Rabbit polyclonal antibody to Shank 2
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human SHANK2. AA range:331-380
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	135
Expiration Date	12 months from date of receipt.



Western Blot analysis of RAW264.7 cells ...



Western blot analysis of lysates from RA...



Western blot analysis of the lysates fro...