



Product Datasheet

TPX2 antibody (orb768155)


Description	Rabbit polyclonal antibody to TPX2
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human DIL-2. AA range:301-350
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	86
Expiration Date	12 months from date of receipt.



Western Blot analysis of Jurkat cells us...



Immunohistochemistry analysis of paraffi...



Western blot analysis of lysates from Ju...