


Product Datasheet

SUN1 antibody (orb76818)

Description	Rabbit monoclonal antibody to SUN1
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ICC, IP, WB
Immunogen	GST-Sun 1 peptide
Target	SUN1
Preservatives	Serum containing 0.1% sodium azide
Note	For research use only
Clonality	Monoclonal
Uniprot ID	O94901
Entrez	23353
Dilution Range	Optimal antibody dilution should be determined by titration, however as a guideline try between 1:500-1:2000
Expiration Date	12 months from date of receipt.



Confocal
immunofluorescence
analysis of ...