



Product Datasheet

QSK antibody (orb768211)

Description	Rabbit polyclonal antibody to QSK
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human QSK. AA range:521-570
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	140
Expiration Date	12 months from date of receipt.



Western Blot analysis of LOVO cells using QSK antibody.



Immunohistochemistry analysis of paraffin-embedded tissue sections using QSK antibody.