


Product Datasheet

TAP antibody (orb76823)

Description	Mouse monoclonal antibody to TAP
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IP, WB
Immunogen	TAP1
Target	TAP1
Preservatives	Purified antibody (from supernatant) containing PBS + 0.1% sodium az
Concentration	1mg/ml
Note	For research use only
Application notes	SP2/0
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	Q03518
Entrez	6890
Dilution Range	WB: 2ug/ml staining in MCF7 cell nuclear lysate. Detects a band of approximately 81 kDa (predicted molecular weight: 87 kDa). ICC/IF: 1:100 on paraformaldehyde fixed MCF7 cells, permeabilized with 0.15% Triton
Expiration Date	12 months from date of receipt.



Western blot analysis of using TAP antib...