


Product Datasheet

POTE-14/22 antibody (orb768257)


Description	Rabbit polyclonal antibody to POTE-14/22
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human A26C2/3. AA range:178-227
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	61
Expiration Date	12 months from date of receipt.



Western Blot analysis of AD293 using POT...



Western blot analysis of lysates from MC...



Immunohistochemical analysis of paraffin...