

Product Datasheet

5-LO (phospho-Ser523) antibody (orb768265)

Description	Rabbit polyclonal antibody to 5-LO (phospho-Ser523)
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Synthesized phospho-peptide around the phosphorylation site of human 5-LO (phospho Ser523)
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	77
Expiration Date	12 months from date of receipt.