


Product Datasheet

Transportin1 antibody (orb76827)

Description	Mouse monoclonal antibody to Transportin1
Species/Host	Mouse
Reactivity	Human, Monkey, Mouse, Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, FC, ICC, IHC-P, IP, WB
Immunogen	Recombinant his tagged transportin 1 (Human)
Target	Transportin 1
Preservatives	Purified antibody (from supernatant) containing PBS 0.1% sodium azide
Concentration	1mg/ml
Note	For research use only
Application notes	SP2/0
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	Q92973 , Q8BFY9
Entrez	238799 , 3842
Dilution Range	Optimal antibody dilution should be determined by titration, however as a guideline try;Flow Cyt: Use 1µg for 106 cells.ICC/IF: Use a concentration of 5 µg/mlIHC-P: Use a concentration of 1 µg/ml. Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol.WB: 1:8000. Detects a band of approximately 93 kDa (predicted molecular weight: 107 kDa)
Expiration Date	12 months from date of receipt.



Western blot
analysis of
Hela lysate
usi...