

Product Datasheet

TRF1 antibody (orb76829)

Description	Mouse monoclonal antibody to TRF1
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, ICC, WB
Immunogen	Human TRF1
Target	TRF1
Preservatives	Purified antibody containing PBS 0.1% sodium azide
Concentration	0.5mg/ml
Note	For research use only
Application notes	Sp2/0 Ag.14
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	P54274
Entrez	7013
Dilution Range	Optimal antibody dilution should be determined by titration, however as a guideline use at 1:100-1:200 for western blot IF. Use Hela cells as positive control
Expiration Date	12 months from date of receipt.