



Product Datasheet

Fyn (phospho-Tyr530) antibody (orb768290)


Description	Rabbit polyclonal antibody to Fyn (phospho-Tyr530)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human Fyn around the phosphorylation site of Tyr530. AA range:488-537
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	60
Expiration Date	12 months from date of receipt.



Western Blot analysis of LOVO + UV cells...



Western blot analysis of 3T3 lysis using...



Immunohistochemistry analysis of paraffi...