


Product Datasheet

Y14 antibody (orb76831)

Description	Mouse monoclonal antibody to Y14
Species/Host	Mouse
Reactivity	Human, Xenopus
Conjugation	Unconjugated
Tested Applications	ELISA, FC, ICC, IHC-P, IP, WB
Immunogen	Full length native protein (purified) (Human)
Target	Y14
Preservatives	Purified antibody (from supernatant) containing PBS 0.1% sodium azide
Concentration	1mg/ml
Note	For research use only
Application notes	Balb/c Sp2/0 cell line
Isotype	IgG2b
Clonality	Monoclonal
Uniprot ID	Q9BRP8
Entrez	9939
Dilution Range	Optimal antibody dilution should be determined by titration, however as a guideline try; Flow Cyt: Use 1-2µg for 10 ⁶ cells IHC-P: Use a concentration of 1 µg/ml. Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol. ICC/IF: 1/15. PubMed: 17103222 WB: Use at an assay dependent dilution. Detects a band of approximately 21 kDa (predicted molecular weight: 21 kDa). PubMed: 17103222
Expiration Date	12 months from date of receipt.



Confocal
immunofluorescence
analysis of ...