



Product Datasheet

Puratrophin 1 antibody (orb768321)


Description	Rabbit polyclonal antibody to Puratrophin 1
Species/Host	Rabbit
Reactivity	Human, Monkey
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human PLEKHG4. AA range:654-703
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	135
Expiration Date	12 months from date of receipt.



Western blot analysis of the lysates from PLEKHG4 cells.



Western blot analysis of lysates from COS cells transfected with PLEKHG4.



Immunohistochemical analysis of paraffin-embedded tissue sections.