



Product Datasheet

Mnk2 antibody (orb768500)

Description	Rabbit polyclonal antibody to Mnk2
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Synthesized peptide derived from the Internal region of human Mnk2.
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	51
Expiration Date	12 months from date of receipt.



Western Blot analysis of HeLa cells using Mnk2 antibody.



Immunohistochemical analysis of paraffin-embedded tissue sections using Mnk2 antibody.