



Product Datasheet

GRB10 (phospho-Tyr67) antibody (orb768506)


Description	Rabbit polyclonal antibody to GRB10 (phospho-Tyr67)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human GRB10 around the phosphorylation site of Tyr67. AA range:33-82
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	67
Expiration Date	12 months from date of receipt.



Western Blot analysis of COLO cells using insulin.



Immunohistochemical analysis of paraffin-embedded tissue sections.



Enzyme-Linked Immunosorbent Assay (Phospho-ELISA)