



Product Datasheet

NMDA zeta 1 antibody (orb768526)


Description	Rabbit polyclonal antibody to NMDA zeta 1
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human NMDAR1. AA range:856-905
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	130
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of NIH/3T3 c...



Immunohistochemistry analysis of paraffi...



Western blot analysis of lysates from Ju...